

# Ana-Chek® CREATINE KINASE - NAC

METHOD – UV KINETIC  
PRODUCT CODE – LC05




## INSTRUCTIONS FOR USE

**INTENDED USE:** Test for estimation of CK - NAC in serum / plasma using UV Kinetic method.

### SUMMARY AND PRINCIPLE

Ana-Chek® CK-NAC is a reagent set for determination of Creatine Kinase activity in human serum/plasma based on UV kinetic method. Ana-Chek® CK-NAC is a two liquid reagent system.



### KIT COMPONENTS

Reagent 1 : CK-NAC Reagent 1  
Reagent 2 : CK-NAC Reagent 2

### REAGENT PREPARATION, STORAGE & STABILITY

Reagent R1 and R2 are ready to use liquid reagents. Mix the reagent R1 and R2 in the ratio of 4:1 respectively to prepare the desired volume of working reagent prior to use. Do not shake vigorously. The working reagent is stable for 14 days at 2-8 °C. The reagent kit should be stored at 2 - 8 °C and is stable till the expiry date indicated on the label.

### PRECAUTIONS & HANDLING

The reagents/samples should be handled by qualified personnel only. Discard reagent/sample as per good laboratory practices and local regulatory requirements. Read the instructions given on the labels and instructions for use carefully before using the kit. The kit is intended for in-vitro diagnostic use only. Don't freeze the reagent. Do not shake the reagent vigorously. Discard the reagent if the absorbance of the reagent exceeds 0.700 O.D. against D/W at 340 nm. Contamination of the reagent should be avoided.

### TEST PARAMETERS

|                          |                |                        |         |
|--------------------------|----------------|------------------------|---------|
| Name                     | CK-NAC         | Reagent Volume         | 1000 µl |
| Reaction Type            | UV Kinetic (↑) | Sample Volume          | 40 µl   |
| Wavelength Primary       | 340 nm         | Incubation Temperature | 37 °C   |
| Flow Cell Temp.          | 37 °C          | Delay Time             | 180 sec |
| Instrument Blank setting | D.W.           | Read Time              | 120 sec |
| Reagent Blank Abs Limit  | <0.700         | Factor                 | 4127    |
| Linearity                | 2000 IU/L      |                        |         |

### MATERIALS REQUIRED BUT NOT PROVIDED

Test tubes, Micropipette with tips, Analyzer, Controls, Incubation chamber.

### SPECIMEN COLLECTION & PRESERVATION

Blood should be collected in a clean dry container. Although serum is preferred, plasma with heparin or EDTA can be used. The assay should

be carried on the same day as far as possible. Serum or plasma are stable for 1 week at 4 °C and 1 month at -10 °C. Avoid use of haemolysed and grossly contaminated samples.

### COMPONENTS OF REAGENT

| Component                                                       | Concentration |
|-----------------------------------------------------------------|---------------|
| Imidazole Buffer, pH 6.7                                        | 100 mmol/l    |
| NAC                                                             | 20 mmol/l     |
| Glucose                                                         | 20 mmol/l     |
| Magnesium Acetate                                               | 10 mmol/l     |
| Hexokinase                                                      | >4500 IU/L    |
| ADP                                                             | 2 mmol/l      |
| AMP                                                             | 5 mmol/l      |
| Di (adenosine -5) - pentaphosphate                              | 10 mol/l      |
| NADP                                                            | 2 mmol/l      |
| EDTA                                                            | 2 mmol/l      |
| Stabilizers and inactive ingredients and surface-active agents. | -             |

### ASSAY PROCEDURE

|                                                                                                                                                          |         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------|---------|
|                                                                                                                                                          | Test    |
| Reagent                                                                                                                                                  | 1000 µl |
| Serum / Plasma                                                                                                                                           | 40 µl   |
| Mix the reagent and sample in the above-mentioned ratio and start the stop watch.                                                                        |         |
| Aspirate reaction mixture into flow cell and record the absorbance at 180 <sup>th</sup> , 210 <sup>th</sup> , 240 <sup>th</sup> , 270 <sup>th</sup> sec. |         |

### CALCULATION

Absorbance/min = Absorbance per 30 sec. x 2  
CK – NAC Activity (IU/L) = Δ Absorbance / min x 4127

### REFERENCE VALUES FOR NORMAL PEOPLE

Men : <190 IU/L at 37 °C.  
Women : <165 IU/L at 37 °C.  
Babies : <325 IU/L at 37 °C.  
Children : <225 IU/L at 37 °C.

### PERFORMANCE CHARACTERISTICS

**Measuring Range:** The assay is linear between 10 – 2000 IU/L. If the CK-NAC value exceeds linearity limit (2000 IU/L), dilute the specimen suitably with normal saline and repeat the assay. In that case, assay value should be multiplied with the dilution factor to obtain correct CK-NAC value of the specimen.

**Interference:** There is no significant interference in samples containing upto 20 mg/dl of bilirubin and 200 mg/dl of haemoglobin.

**Precision:** Precision studies has been carried out using quality control sera as shown below:

| (n=10)            | Within Run  |           |      | Between Run |           |      |
|-------------------|-------------|-----------|------|-------------|-----------|------|
| Specimen Material | Mean (IU/L) | SD (IU/L) | CV % | Mean (IU/L) | SD (IU/L) | CV % |
| Low Value Serum   | 141.47      | 1.72      | 1.2  | 147.03      | 1.96      | 1.3  |
| High Value Serum  | 417.0       | 1.45      | 0.3  | 427.7       | 2.63      | 0.6  |

Note: We recommend all the laboratories to establish its own accuracy and precision data.

## QUALITY CONTROL

Inclusion of a normal value and abnormal value chemistry control serum in each test run ensures optimum quality control. Consistent use of same type and methodology of control serum provides between run precision and accuracy data for CK-NAC. We recommend to produce such data on daily basis for greater accuracy in assay system which include reagents, instrument, apparatus and operator.

## PRECAUTIONS

1. The working reagent is considered unsatisfactory and should not be used if its absorbance exceeds 0.700 at 340 nm against distilled water.
2. If CK-NAC value exceeds 2000 IU/L then dilute the specimen suitably with normal saline & repeat the assay.
3. Avoid using haemolysed serum since red blood cells may release enzymes and intermediates such as ATP and G6P which may interfere.
4. CK-NAC Activity in serum may be elevated in patients receiving intramuscular injection upto one week prior to sample collection. Strenuous or unusual physical exercise may also cause elevated CK-NAC activity.

## BIBLIOGRAPHY

1. Stein W. Laboratory Diagnosis of Acute myocardial infarction. Darmstadt: GIT verlag, 1988:34-37.
2. Horder M, Elser RC, Gerhardt M *et al*, Approved Recommendation on IFCC Methods for the Measurement of catalytic Concentration of Enzymes. part7. IFCC method for Creatine Kinase. Eur. J. Clin. Chem. Clin. Biochem., 1991; 29:435-456.
3. Tietz NW, ed. Clinical Guide to laboratory Tests, 3rd ed. Philadelphia, Pa: W.B. Saunders Company., 1995:181-181.
4. Scandinavian society for Clinical chemistry and clinical Physiology scand J.Clin.Lab.Invest.,37,711(1976)
5. Meattini F, Giannini G. and Tarli P, clin.chem; 24,3(1978)

| Symbol                                                                              | Explanation           | Symbol                                                                                | Explanation                  |
|-------------------------------------------------------------------------------------|-----------------------|---------------------------------------------------------------------------------------|------------------------------|
|  | Manufactured By       |  | In Vitro Diagnostic Use      |
|  | Lot Number            |  | Read Instructions Before Use |
|  | Catalogue Number      |  | Storage Temperature          |
|  | Manufacturing Date    |  | Number of Tests / Volume     |
|  | Expiry Date           |  | Do Not Reuse                 |
|  | Protect from Sunlight |  | Keep Dry                     |

**ANAMOL LABORATORIES PVT. LTD.**  
61, Genesis Industrial Township, Kolgaon,  
Palghar – 401 404. India.  
Customer Care & WhatsApp: +91-9823388695.

admin@anamollabs.com  
exports@anamollabs.com  
www.anamollabs.com

ISO 9001: 2015  
ISO 13485: 2016  
GMP  
CE